FLEXIDY®

ISTRUZIONI D'USO FLEXIDY®

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Trade Name:	FLEXIDY®
Chemical Composition:	Ethylene-Vinyl-Acetate (EVA) Copolymer.
Molecular structure :	Low Crystallinity Thermoplastic Material.
Class according to EEC/93/42: I.	
Product form:	Cylinders.
Use indications:	Positioners and Mouth-guards for sporting.

FLEXIDY® is a Thermoplastic Technopolymer to be worked by MG-NEWPRESS® (Thermopress) or FLEXIPRESS® (Manual Injector) heat-molding systems. The indications set forth herein illustrate the correct operational parameters to follow to guarantee an optimal result according to the system used.

Guide to Pressing by MG-NEWPRESS® Thermo-Injection process

FLEXIDY® Type	"80" (Hard)	"65" (Medium)
Temperature of the Oven	160 °C	120 °C
Permanence of material in the Oven	15 minutes	15 minutes
Temperature of Mould/Muffle	23 °C about.	23 °C about.
Injection Pressure	3 Atm (Bar)	3 Atm (Bar)
Post-Injection Pressure: (2 seconds after Injection)	6 Atm (Bar)	6 Atm (Bar)
Initial Cooling: Permanence with Piston and Oven On	10 minutes	10 minutes
Final Cooling: Permanence of Muffle in cool water	25 minutes	25 minutes

Guide to Pressing by FLEXIPRESS® Manual System

FLEXIDY® Type	"80" (Hard)	"65" (Medium)
Temperature of water	100 °C	100 °C
Permanence of material in water	45 minutes	5 minutes
Temperature of mould/Muffle	100 °C	40 °C
Material plastification method	Immersion of Piston and Muffle	Immersion of Piston only
Post-Injection pressure (screw screwing after injection)	YES	NO
Final Cooling: permanence of Muffle in cool water	20 minutes	30 minutes

IMPORTANTE:

for detailed information on working procedures, consult MG-NEWPRESS® User's Manual, sections "General principles on MG-NEWPRESS® injection system" (pag.7) and "Use with FLEXIDY®" (pag.9), or FLEXIPRESS® User's Manual.

WARNINGS:

1. Before plasticizing procedure carefully isolate Muffle by means of a common plaster/resin isolator.
2. Plasticize FLEXIDY® only with MG-NEWPRESS® or with FLEXIPRESS®.
3. Do not use higher temperatures than the recommended ones as there is a risk that the technopolymer may suffer heat decomposition.
4. Always and only utilize virgin material, do not mix with leftover pieces from previous plasticizations which may pollute the material, compromising its physical and mechanical characteristics
5. Do not use hot water to wash FLEXIDY®, only cold water must be used.

STORAGE:

Store in a dry and airy place at room temperature.

CLEANING:

Products made out of FLEXIDY® may be cleaned and disinfected using common specific cleansers.

CE

QUATTROTI

QuattroTi S.r.l. - Via Vittorio Veneto, 34
22070 Rovello Porro (CO) Italia
Tel. +39 02.96750972 - Fax +39 02.96751450